



Through Vial Impedance Spectroscopy (TVIS)

A new PAT for freeze-drying cycle development

Smi Pharmaceutical Freeze Drying Technology

13th - 14th June 2018

Holiday Inn Kensington Forum, London, United Kingdom

Prof. Geoff Smith


Leicester School of Pharmacy, De Montfort University, United Kingdom

Thanks to **GEA Pharma Systems** for financial support to attend the conference


Outline



- Introduction
 - Description of Measurement System
 - Theory (Dielectric loss mechanisms)
- TVIS Applications
 - Ice Formation and phase separation in freezing
 - Temperature calibration and prediction
 - Drying rate estimation
 - Heat transfer coefficient calculation
 - End-point determination
 - Dry Layer resistance and collapse (time permitting)
- Acknowledgements



Yowwares Jeeraruangrattana
GPO Thailand



TVIS pass through on
GEA Lyophil dryer, Hurth, Cologne

Introduction to the TVIS System






- Impedance spectroscopy characterizes the ability of materials to conduct electricity under an applied an oscillating voltage (of varying frequency)
- Impedance measurements **across a vial** rather than **within the vial**
- Hence **“Through Vial Impedance Spectroscopy”**
- Features
 - Single vial “non-product invasive”
 - Both freezing and drying characterised in a single technique
 - Non-perturbing to the packing of vials
 - Stopper mechanism unaffected




Through Vial Impedance Spectroscopy (TVIS)

Description of Measurement System


Freeze drying chamber



LyoView™ analysis software




LyoDEA™ measurement software



**TVIS system
(I to V convertor)**




TVIS Measurement System



- The designs of various vials that have been modified with copper foil electrodes (10 mm in height and 3 mm from the base of each container)
 - 20 mm crimp-neck vial with 10 ml nominal capacity
 - 20 mm crimp-neck vial with 5 ml nominal capacity
 - screw-neck vial with 5 ml nominal capacity

- The different styles of a bespoke pass-through for TVIS systems



- Connected via the manifold hose on the outside of the dryer
- Connected via the port on top left side of the door on the dryer
- Connected to a port on the top of the drying chamber




Through Vial Impedance Spectroscopy (TVIS)

Dielectric Loss Mechanisms


Measurement vial



Electronic polarization
distortion of electrons
relative to the nuclei




Atomic polarization
distortion of nuclei across
a heteroatom bond by
stretching and bending




Instantaneous polarization dominant
mechanism in the glass wall

Electronic polarization
distortion of electrons
relative to the nuclei



Atomic polarization
distortion of nuclei across
a heteroatom bond by
stretching and bending





Instantaneous polarization dominant
mechanism in the glass wall


TVIS response for empty vial

TVIS response for empty vial


Measurement vial




Electronic polarization
distortion of electrons
relative to the nuclei



Atomic polarization
distortion of nuclei across
a heteroatom bond by
stretching and bending



**Instantaneous polarization dominant
mechanism in the glass wall**




**Space charge polarization
(weak frequency dependence)**


TVIS response for empty vial

MW (space-charge) polarization at glass wall – sample interface


Electronic polarization
distortion of electrons
relative to the nuclei



Atomic polarization
distortion of nuclei across
a heteroatom bond by
stretching and bending





**Instantaneous polarization dominant
mechanism in the glass wall**




**Space charge polarization
(weak frequency dependence)**

TVIS response for empty vial


Measurement vial




Electronic polarization
distortion of electrons
relative to the nuclei



Atomic polarization
distortion of nuclei across
a heteroatom bond by
stretching and bending




**Instantaneous polarization dominant
mechanism in the glass wall**



**Space charge polarization
(weak frequency dependence)**

**MW (space-charge) polarization at glass
wall – sample interface**

TVIS response for empty vial




Proton-hopping

Conduction of protons in liquid water occurs
through the Grotthuss "hop-turn" mechanism

Conductivity in pure water



Dipolar polarization

re-orientation/alignment of
permanent dipoles in liquid
water (Debye-like
relaxation)



**Polarization mechanisms in liquid water
(relaxation time, $\tau \sim 9$ ps at 20°C)**

TVIS response for liquid water




Electronic polarization
distortion of electrons
relative to the nuclei



Atomic polarization
distortion of nuclei across
a heteroatom bond by
stretching and bending

**Instantaneous polarization dominant
mechanism in the glass wall**





**Space charge polarization
(weak frequency dependence)**


**MW (space-charge) polarization at glass
wall – sample interface**

TVIS response for empty vial

Measurement vial




Electronic polarization
distortion of electrons
relative to the nuclei



Atomic polarization
distortion of nuclei across
a heteroatom bond by
stretching and bending



**Instantaneous polarization dominant
mechanism in the glass wall**




**Space charge polarization
(weak frequency dependence)**


**MW (space-charge) polarization at glass
wall – sample interface**

TVIS response for empty vial

Dominant at $T > 235$ K (approx. -40 °C)
Generation/migration of L- and D- orientation defects
in ice Ih




Dominant at $T < 235$ K (approx. -40 °C)
Generation/migration of H_3O^+ / OH^- ion pairs (ionic
defects) in ice Ih (similar to the Grotthus mechanism)




Polarization mechanism in ice

TVIS response for frozen water (ice)


Electronic polarization
distortion of electrons
relative to the nuclei



Atomic polarization
distortion of nuclei across
a heteroatom bond by
stretching and bending



**Instantaneous polarization dominant
mechanism in the glass wall**





**Space charge polarization
(weak frequency dependence)**

**MW (space-charge) polarization at glass
wall – sample interface**

TVIS response for empty vial

Frozen Water and Dielectric Relaxation of Ice



- I. : The polarization of the water dipole in liquid water at 20 °C, with a dielectric loss peak frequency of ~ 18 GHz



Frozen Water and Dielectric Relaxation of Ice





- I. : The polarization of the water dipole in liquid water at 20 °C, with a dielectric loss peak frequency of ~ 18 GHz
- II. : The Maxwell-Wagner (MW) polarization of the glass wall of the TVIS vial at 20 °C, with a dielectric loss peak frequency of 17.8 kHz



Frozen Water and Dielectric Relaxation of Ice





- I. : The polarization of the water dipole in liquid water at 20 °C, with a dielectric loss peak frequency of ~ 18 GHz
- II. : The Maxwell-Wagner (MW) polarization of the glass wall of the TVIS vial at 20 °C, with a dielectric loss peak frequency of 17.8 kHz
- III. : The dielectric polarization of ice at -20 °C, with a dielectric loss peak frequencies of 2.57 kHz



Frozen Water and Dielectric Relaxation of Ice



- I. : The polarization of the water dipole in liquid water at 20 °C, with a dielectric loss peak frequency of ~ 18 GHz
- II. : The Maxwell-Wagner (MW) polarization of the glass wall of the TVIS vial at 20 °C, with a dielectric loss peak frequency of 17.8 kHz
- III. : The dielectric polarization of ice at -20 °C, with a dielectric loss peak frequencies of 2.57 kHz
- IV. : The dielectric polarization of ice at -40 °C with a dielectric loss peak frequencies of 537 Hz.




Through Vial Impedance Spectroscopy (TVIS)

Applications

Through Vial Impedance Spectroscopy (TVIS)



Monitoring **Phase Behaviour**
(ice nucleation temperature
and solidification end points
by using F_{PEAK})





Through Vial Impedance Spectroscopy (TVIS)



Monitoring **Phase Behaviour**
(ice nucleation temperature
and solidification end points
by using F_{PEAK})


F_{PEAK} temperature calibration
for **predicting temperature** of
the product in primary drying




Through Vial Impedance Spectroscopy (TVIS)




Monitoring **Phase Behaviour**
(ice nucleation temperature
and solidification end points
by using F_{PEAK})



F_{PEAK} temperature calibration
for **predicting temperature** of
the product in primary drying




Surrogate **drying rate**
(from $\frac{dC''_{PEAK}}{dt}$)




Through Vial Impedance Spectroscopy (TVIS)




Monitoring **Phase Behaviour**
(ice nucleation temperature
and solidification end points
by using F_{PEAK})



F_{PEAK} temperature calibration
for **predicting temperature** of
the product in primary drying



Surrogate **drying rate**
(from $\frac{dC''_{PEAK}}{dt}$)



C' (~ 100 kHz) is highly sensitive to low ice volumes; therefore it could be used for determination **end point** of primary drying



TVIS Response Surface (3D-Plot)



Imaginary Part of Capacitance


Real Part of Capacitance

Annealing = Re-heating and Re-cooling






Phase Separation in Freezing Step

Water (frozen)




5%w/v Lactose solution (frozen)



Dielectric loss spectrum


- Data analysing software (*LyoView™*) identifies the peak frequency (F_{PEAK}) and peak amplitude (C''_{PEAK}) in the imaginary part of the capacitance spectrum



Through Vial Impedance Spectroscopy (TVIS)



Temperature calibration & Temperature prediction

Temperature Calibration





Polynomial coefficient from $\text{Log } F_{PEAK}$ – temperature calibration

	a	b	c
TE	-1.02	30.1	-114
BE	-5.43×10^{-1}	26.7	-109





Temperature Prediction





Through Vial Impedance Spectroscopy (TVIS)

Drying rate estimation


Drying Rate Estimation



Drying Rate Estimation




Drying Rate Estimation



- Drying rate during the steady state

$$\text{Drying rate } \left(\frac{\Delta m}{\Delta t} \right) = \rho_i \cdot A \cdot \frac{h_{(t_1)} - h_{(t_2)}}{t_2 - t_1}$$

Drying Rate Estimation




- Drying rate during the steady state

$$\text{Drying rate } \left(\frac{\Delta m}{\Delta t} \right) = \rho_i \cdot A \cdot \frac{h_{(t_1)} - h_{(t_2)}}{t_2 - t_1}$$

Ice density (ρ_i) at -32°C = $0.920 \text{ g}\cdot\text{cm}^{-3}$
 (Calculated ice temperature between T_i & T_b)

Drying Rate Estimation



- Drying rate during the steady state


$$\text{Drying rate } \left(\frac{\Delta m}{\Delta t} \right) = \rho_i \cdot A \cdot \frac{h_{(t_1)} - h_{(t_2)}}{t_2 - t_1}$$

Ice density (ρ_i) at -32°C = $0.920 \text{ g} \cdot \text{cm}^{-3}$

(Calculated ice temperature between T_i & T_b)

Internal vial diameter (VC010-20C) = 2.21 cm

Drying Rate Estimation



- Drying rate during the steady state

$$\text{Drying rate } \left(\frac{\Delta m}{\Delta t} \right) = \rho_i \cdot A \cdot \frac{h_{(t1)} - h_{(t2)}}{t_2 - t_1}$$


Ice density (ρ_i) at -32°C = $0.920 \text{ g}\cdot\text{cm}^{-3}$

(Calculated ice temperature between T_i & T_b)

Internal vial diameter (VC010-20C) = 2.21 cm

Cross-section area (A) = 3.80 cm^2

Drying Rate Estimation



- Drying rate during the steady state

$$\text{Drying rate } \left(\frac{\Delta m}{\Delta t} \right) = \rho_i \cdot A \cdot \frac{h_{(t1)} - h_{(t2)}}{t_2 - t_1}$$

Ice density (ρ_i) at -32°C = $0.920 \text{ g}\cdot\text{cm}^{-3}$


(Calculated ice temperature between T_i & T_b)

Internal vial diameter (VC010-20C) = 2.21 cm

Cross-section area (A) = 3.80 cm^2

Ice height at 2 h ($h_{(2h)}$) = 19.94 mm

Drying Rate Estimation



- Drying rate during the steady state

$$\text{Drying rate } \left(\frac{\Delta m}{\Delta t} \right) = \rho_i \cdot A \cdot \frac{h_{(t_1)} - h_{(t_2)}}{t_2 - t_1}$$

Ice density (ρ_i) at -32°C = $0.920\text{ g}\cdot\text{cm}^{-3}$

(Calculated ice temperature between T_i & T_b)


Internal vial diameter (VC010-20C) = 2.21 cm

Cross-section area (A) = 3.80 cm^2

Ice height at 2 h ($h_{(2\text{ h})}$) = 19.94 mm

Ice height at 2.8 h ($h_{(2.8\text{ h})}$) = 18.98 mm

Drying Rate Estimation



- Drying rate during the steady state

$$\text{Drying rate } \left(\frac{\Delta m}{\Delta t} \right) = \rho_i \cdot A \cdot \frac{h_{(t_1)} - h_{(t_2)}}{t_2 - t_1}$$

Ice density (ρ_i) at -32°C = $0.920 \text{ g} \cdot \text{cm}^{-3}$

(Calculated ice temperature between T_i & T_b)

Internal vial diameter (VC010-20C) = 2.21 cm


Cross-section area (A) = 3.80 cm^2

Ice height at 2 h ($h_{(2h)}$) = 19.94 mm

Ice height at 2.8 h ($h_{(2.8h)}$) = 18.98 mm

$$\begin{aligned} \text{Drying rate} &= 0.920 \text{ g} \cdot \text{cm}^{-3} \times 3.80 \text{ cm}^2 \times \frac{(19.94 - 18.98) \times 10^{-1} \text{ cm}}{(2.8 - 2.0) \text{ h}} \\ &= \mathbf{0.42 \text{ g} \cdot \text{h}^{-1}} \end{aligned}$$

Drying Rate Estimation



- Drying rate during the steady state

$$\text{Drying rate } \left(\frac{\Delta m}{\Delta t} \right) = \rho_i \cdot A \cdot \frac{h_{(t_1)} - h_{(t_2)}}{t_2 - t_1}$$

Ice density (ρ_i) at -32°C = $0.920 \text{ g} \cdot \text{cm}^{-3}$

(Calculated ice temperature between T_i & T_b)

Internal vial diameter (VC010-20C) = 2.21 cm

Cross-section area (A) = 3.80 cm^2

Ice height at 2 h ($h_{(2h)}$) = 19.94 mm

Ice height at 2.8 h ($h_{(2.8h)}$) = 18.98 mm

TVIS parameters used for determination:

$$\frac{\Delta m}{\Delta t} = 0.42 \text{ g} \cdot \text{h}^{-1}$$

$$T_b = -29.8^{\circ}\text{C}$$

$$\begin{aligned} \text{Drying rate} &= 0.920 \text{ g} \cdot \text{cm}^{-3} \times 3.80 \text{ cm}^2 \times \frac{(19.94 - 18.98) \times 10^{-1} \text{ cm}}{(2.8 - 2.0) \text{ h}} \\ &= \mathbf{0.42 \text{ g} \cdot \text{h}^{-1}} \end{aligned}$$

Through Vial Impedance Spectroscopy (TVIS)

Heat Transfer Coefficient Determination

Heat Transfer Coefficient (K_v) Determination



Parameters	TVIS
Drying rate at steady state (g/h) (2-2.8 h into primary drying)	0.42
Shelf Temperature, T_s (K)	273.3
Vial's base Temperature, T_b (K)	243.3

Heat Transfer Coefficient (K_v) Determination



Parameters	TVIS
Drying rate at steady state (g/h) (2-2.8 h into primary drying)	0.42
Shelf Temperature, T_s (K)	273.3
Vial's base Temperature, T_b (K)	243.3

$$L \frac{\Delta m}{\Delta t} = A_e K_v (T_s - T_b) \Rightarrow K_v = \frac{L \frac{\Delta m}{\Delta t}}{A_e (T_s - T_b)}$$

L is the latent heat of sublimation of ice (2844 J·g⁻¹ or 679.7 cal·g⁻¹) and A_e is external cross-sectional area of the base of the TVIS vial (4.62 cm²)

Heat Transfer Coefficient (K_v) Determination



Parameters	TVIS
Drying rate at steady state (g/h) (2-2.8 h into primary drying)	0.42
Shelf Temperature, T_s (K)	273.3
Vial's base Temperature, T_b (K)	243.3

$$L \frac{\Delta m}{\Delta t} = A_e K_v (T_s - T_b) \Rightarrow K_v = \frac{L \frac{\Delta m}{\Delta t}}{A_e (T_s - T_b)}$$

L is the latent heat of sublimation of ice (2844 J·g⁻¹ or 679.7 cal·g⁻¹) and A_e is external cross-sectional area of the base of the TVIS vial (4.62 cm²)

$$K_v(270 \mu\text{bar}) = \frac{L \frac{\Delta m}{\Delta t}}{A_e (T_s - T_b)}$$


$$\begin{aligned} &= \frac{679.7 \text{ cal} \cdot \text{g}^{-1} \times 0.42 \text{ g} \cdot \text{h}^{-1}}{4.62 \text{ cm}^2 \times (273.3 - 243.3) \text{ K}} \\ &= 2.06 \text{ cal} \cdot \text{h}^{-1} \cdot \text{cm}^{-2} \cdot \text{K}^{-1} \\ &= 5.73 \times 10^{-4} \text{ cal} \cdot \text{s}^{-1} \cdot \text{cm}^{-2} \cdot \text{K}^{-1} \end{aligned}$$

$$K_v(270 \mu\text{bar}) = 5.73 \times 10^{-4} \text{ cal} \cdot \text{s}^{-1} \cdot \text{cm}^{-2} \cdot \text{K}^{-1}$$

Heat Transfer Coefficient (K_v) Determination



Parameters	TVIS
Drying rate at steady state (g/h) (2-2.8 h into primary drying)	0.42
Shelf Temperature, T_s (K)	273.3
Vial's base Temperature, T_b (K)	243.3



$$L \frac{\Delta m}{\Delta t} = A_e K_v (T_s - T_b) \Rightarrow$$

$$K_v = \frac{L \frac{\Delta m}{\Delta t}}{A_e (T_s - T_b)}$$

$$K_v(270 \mu\text{bar}) = \frac{L \frac{\Delta m}{\Delta t}}{A_e (T_s - T_b)}$$

L is the latent heat of sublimation of ice ($2844 \text{ J} \cdot \text{g}^{-1}$ or $679.7 \text{ cal} \cdot \text{g}^{-1}$) and A_e is external cross-sectional area of the base of the TVIS vial (4.62 cm^2)





$$\begin{aligned} &= \frac{679.7 \text{ cal} \cdot \text{g}^{-1} \times 0.42 \text{ g} \cdot \text{h}^{-1}}{4.62 \text{ cm}^2 \times (273.3 - 243.3) \text{ K}} \\ &= 2.06 \text{ cal} \cdot \text{h}^{-1} \cdot \text{cm}^{-2} \cdot \text{K}^{-1} \\ &= 5.73 \times 10^{-4} \text{ cal} \cdot \text{s}^{-1} \cdot \text{cm}^{-2} \cdot \text{K}^{-1} \end{aligned}$$

$$K_v(270 \mu\text{bar}) = 5.73 \times 10^{-4} \text{ cal} \cdot \text{s}^{-1} \cdot \text{cm}^{-2} \cdot \text{K}^{-1}$$





Through Vial Impedance Spectroscopy (TVIS)

End-point Determination





End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)






End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)






End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)







End-Point Determination (T_s -15 °C and P_c 400 μ bar)






End-Point Determination (T_s -15 °C and P_c 400 μ bar)



End-Point Determination (T_s -15 °C and P_c 400 μ bar)




End-Point Determination (T_s -15 °C and P_c 400 μ bar)




Summary

- End point determine by high frequency (i.e. 100 kHz for ice) real part capacitance (i.e. 100 kHz for ice)



Summary

- End point determine by high frequency (i.e. 100 kHz for ice) real part capacitance (i.e. 100 kHz for ice)





Through Vial Impedance Spectroscopy (TVIS)


Product Resistance Determination

Product Resistance (R_p) Determination

$$\hat{R}p = \left(\frac{P_{ICE} - P_{CHAMBER}}{dm/dt} \right) \cdot A_p$$







Lactose freeze-dried matrix




Product Resistance (R_p) Determination

$$\hat{R}_p = \left(\frac{P_{ICE} - P_{CHAMBER}}{dm/dt} \right) \cdot A_p$$




Top layer
Fine pores



Middle layer
Micro-collapse

500X






Bottom layer
Full collapse



Through Vial Impedance Spectroscopy (TVIS)

Collapse Phenomena


5% Sucrose




$T_{(TC)} = -37.8^\circ\text{C};$
 $T_{(FPEAK)} = -38.2^\circ\text{C}$




$T_{(TC)} = -38.2^\circ\text{C};$
 $T_{(FPEAK)} = -38.2^\circ\text{C}$




$T_{(TC)} = -34.9^\circ\text{C};$
 $T_{(FPEAK)} = -35.6^\circ\text{C}$





$T_{(TC)} = -35.2^\circ\text{C};$
 $T_{(FPEAK)} = -36.0^\circ\text{C}$





$T_{(TC)} = -33.0^\circ\text{C};$
 $T_{(FPEAK)} = -33.6^\circ\text{C}$




5% Sucrose in 0.26% NaCl




$T_{(TC)} = -39.5$ °C;
 $T_{(F_{PEAK})} = -37.5$ °C




$T_{(TC)} = -40.2$ °C;
 $T_{(F_{PEAK})} = -37.8$ °C




$T_{(TC)} = -38.0$ °C;
 $T_{(F_{PEAK})} = -37.0$ °C






$T_{(TC)} = -38.1$ °C;
 $T_{(F_{PEAK})} = -37.0$ °C



$T_{(TC)} = -36.8$ °C;
 $T_{(F_{PEAK})} = -36.1$ °C



5% Sucrose +0.55% NaCl



Summary results of sugar- salts solutions



Parameters	Unit	S-1	S-2	S-3
Component				
Sucrose	%w/v	5	5	5
Sodium chloride	%w/v	-	0.26	0.55
T'_g (DSC)	°C	-33.7	-38.0	-41.8
Temperature stabilization period (Before ramping shelf temperature)		III to IV	III to IV	III to IV
	h	23.3 - 23.6	23.3 - 23.6	23.1 - 23.6
$T_{(F_{PEAK})}$	°C	-36.0 to -35.6	-37.0	-38.3 to -37.7
Surrogate drying rate ($d\hat{C}_{PEAK}/dt$)	pF/h	0.0380	0.0644	0.0864
Structural state		Non-collapse	micro-collapse?	Micro-collapse
		VII to VIII	VII to VIII	VII to VIII
Temperature stabilization period (After shelf temperature constant)	h	26.0 - 26.8	26.0 – 26.4	26.0 -26.2
	°C	-26.4 to -25.6	-30.1 to -30.0	-32.8 to -32.7
Surrogate drying rate ($d\hat{C}_{PEAK}/dt$)	pF/h	0.2243	0.1227	0.1569
Structural State*		Micro-collapse	Collapse	Collapse

5% Sucrose in Salt Solutions

VI. 23.90 h



$T_{(TC)} = -36.8\text{ }^{\circ}\text{C};$
 $T_{(FPEAK)} = -36.1\text{ }^{\circ}\text{C}$

VIII. 26.40 h



$T_{(TC)} = -31.7\text{ }^{\circ}\text{C};$
 $T_{(FPEAK)} = -30.0\text{ }^{\circ}\text{C}$

VII. 26.03 h




$T_{(TC)} = -31.5\text{ }^{\circ}\text{C};$
 $T_{(FPEAK)} = -30.1\text{ }^{\circ}\text{C}$

IX. 26.67 h



$T_{(TC)} = -31.4\text{ }^{\circ}\text{C};$
 $T_{(FPEAK)} = -27.9\text{ }^{\circ}\text{C}$



5% sucrose in 0.26% NaCl

VI. 24.20 h



$T_{(TC)} = -36.7\text{ }^{\circ}\text{C};$
 $T_{(FPEAK)} = -36.2\text{ }^{\circ}\text{C}$

VIII. 26.20 h



$T_{(TC)} = -34.4\text{ }^{\circ}\text{C};$
 $T_{(FPEAK)} = -32.7\text{ }^{\circ}\text{C}$

VII. 26.03 h




$T_{(TC)} = -34.1\text{ }^{\circ}\text{C};$
 $T_{(FPEAK)} = -32.8\text{ }^{\circ}\text{C}$

IX. 26.90 h



$T_{(TC)} = -34.7\text{ }^{\circ}\text{C};$
 $T_{(FPEAK)} = -32.5\text{ }^{\circ}\text{C}$



5% sucrose in 0.55% NaCl

Limitations

- C''_{PEAK} and F_{PEAK} parameters rely on intimate contact of ice cylinder with glass wall
- $C'(100\text{ kHz})$ parameter does not depend on contact and can be used for end point but relationship between $C'(100\text{ kHz})$ ice constant is non-linear
- Cable length limited to 1m at present
- C-TVIS not compatible with front loading system
- Incompatible with TCs in same TVIS vial (use fibre optic sensors – INFAP)

Future Work

- Development mapping a drying characteristics
 - heat transfer coefficients (K_V)
 - dry layer resistance (R_P)



- TVIS vial – Electrode Coating!
 - Electrodeposition
- Instrument Development
 - Commercial C-TVIS (2018)
 - Non-contact TVIS (2018-19)
 - Micro-well screening
 - Vial clusters in batch FD
 - TVIS - Shuttle (2019-20)

Non-invasive real time information for characterising the freeze drying

Acknowledgements, Recent Projects & Collaborators

- De Montfort University, School of Pharmacy
 - Evgeny Polygalov: co-inventor of TVIS instrument
 - Yowwares Jeeraruangrattana. PhD student
 - Bhaskar Pandya. PhD student
 - Irina Ermolina. Senior Lecturer

